
NARVA

Elektrische Radiometerröhre ERR 2

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Kleine Anodenspannung !

Länge ca. 220 mm

Ø 125 mm

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Tube holder with relaxation generator RHK 4e

The tube holder and relaxation generator from one unit.

When the toggle switch is set to the relaxation generator work can be performed with the cathode ray tube. When the toggle switch is set downwards, the relaxation generator is switched off. The unit is now a tube holder and can be used to display all tubes for teaching purposes from Glimmlampenwerk Cursdorf.

The holder comprises tube socket, two test terminals for ac magnets and two test terminals for cathode ray tube deflection. At the side is the connecting cable for the supply unit.

The colours are:

plate	grid	mass	heating
red	yellow	black	green

Needed to represent the sinusoidal curve of the mains voltage are:

1 tripod (Apolda)

1 BR 2 cathode ray tube (Cursdorf)

1 electromagnet (6.3 volt)

1 power supply unit with the following outputs:

6.3 volts for heating and magnet

200—400 volts plate voltage

50 volts grid voltage

Direct current voltages must be well smoothed.

For experiments with the cathode ray tube a holder can be supplied to prevent the tube base being overloaded and breaking off.

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Support for Tubes R. H. 2 (Frama)

Referring all vakuum-tubes with cathode indirect heated, produced by Glimmlampenwerk of Cursdorf (Thür.)

On the support are installed the hold for the tubes, two holdfasts for alternating current-magnet, two holdfasts to deflect the tubes according to Braun. At the top there are bushes to join the tool for canting these tubes. At the side there is the supplycord for the tool providing the current.

The colours are:

Anode:	Lattice:	Mass:	Heating:
red	yellow	black	green

To represent the sine-curve of the intensity of a current are needed :

1 stand (tripod—Apolda)

1 tube according to Braun

1 tool for canting

1 current-magnet (6,3 Volt)

1 tool for delivery of the current with the following startings.

6,3 Volt for heating and magnet

200—400 Volt anode-tension

—50 Volt lattice-tension

Tensions of continous current are smoothed well.

On experiments with employment of the tubes according to Braun a support can be deliveres to do not overcharge this tube (risk of breaking off.).

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